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
Clinical Case Report / Relato de Caso Clínico

Bionnovation
b i o m e d i c a l

by: **Helder Valiense, DDS, Ms (Ilhéus, Bahia - Brazil)**

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Female Patient, 36 years old, came to the clinic with complains of dissatisfaction with the dental implants. In the radiographic examination revealed loss of alveolar bone structure and Implants threads are exposed and proper osseointegration.

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Treatment proposal: GBR (Guided Bone Regeneration) for both Implants.

Was used Titanium Mesh (Surgitime Titanium), Bovine Bone (Bonefill® Dense) Fine and Medium granulation. The implants functioned as supporting pillars for the titanium mesh, so as to sustain and allowing the adequate and safe accommodation to tent the tissue volume. The insertion of the Surgitime Titanium was to help to protect against the formation of scar tissue at the surgical site. The Bonefill® was placed into the defect in layers: A Medium-grained of Bonefill® was placed just above of the Fine granulation.

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
6 months after graft procedure demonstrates excellent reconstruction of defect and coverage of exposed implant threads with newly formed bone. Reopening the surgical site, removing the titanium mesh and using a provisional restorations on implants.

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